

Plum



The +

- * Available in several packages
- * Sterile eyewash
- * Same composition as eye fluid

PURCHASE PACKAGING

Ref.	Size	Inner	Carton
GANPL4701	200 ml		

DESCRIPTION

- * sterile, 0,9% sodium chloride solution (physiological salt solution, in composition identical with eye liquid)
- * protects cornea and crystalline lens when dust and other foreign materials get into eyes
- * holes on the cup of dust-proof eye-washing head prevent further infection
- * rinsing time: 2 or 5 minutes/bottle (100 ml/minute)

APPLICATIONS

The rinse of eyes if foreign matter gets into eyes

TECHNICAL FEATURES

Features	* 200 ml * rinsing time: 2 or 5 minutes/bottle (100 ml/minute)
Main risk	* complies with medical device regulations protects cornea and crystalline lens when dust and other foreign materials get into eyes

STANDARD(S)

RETAILER STAMP

MAJ 21/06/2022